

# NASO-PHARYNGO-LARYNGOSCOPE

## OP-30

HIGH RESOLUTION IN THE SMALLEST DIAMETER

- Up to 18.000 pixels
- Detailed vision
- High resistance



producing quality

OPTOMIC flexible endoscopes provide excellent optical characteristics for the inspection of nasal passages, the pharynx and the larynx. Both fiberscopes have been designed to facilitate examination efficiency and enhance patient comfort.

# NASO-PHARYNGO-LARYNGOSCOPE **OP-30**



FIBROLUX LED DUO LIGHT SOURCE



FIBROLUX LED LIGHT SOURCE



FIBER OPTIC CABLE



FIBERSCOPE HOLDER



LEAKAGE TEST DEVICE



TWO-HANDED GRIP

	OP-30
Image guide system	18.000 pixels
Outer diameter	3,4mm
Flexible part diameter	3,2mm
Depth of field	2,5-50mm
Deflection angle	130°up/down
Direction of view	0°
Working length	300mm

- High definition images (up to 18.000 pixels).
- Focusing adjustment.
- ETO valve.
- Soakable in disinfectant liquid.
- Light guide connector for Storz and ACMI (optional Wolf).
- Delivered in a protective carrying case.

## OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain

Tel: +34 91 804 44 44

export@optomic.com

www.optomic.com



producing quality

Specifications are subject to change without notice