

# ENT RIGID ENDOSCOPES

## OTO-RHINOSCOPE-LARYNGOSCOPES

OUTSTANDING QUALITY IN ENDOSCOPIC VISION


- High resolution
- Clear images
- High resistance



OPTOMIC offers a complete range of rigid otoscopes, rhinoscope and laryngoscopes for the inspection of the corresponding area. They allow for the diagnosis and serve as a method for minimal invasive surgical intervention. All the ENT OPTOMIC endoscopes are of proven quality and offer an exceptional image brightness to guarantee the satisfaction of the specialist.

## ENT RIGID ENDOSCOPES


# OTO-RHINOSCOPE-LARYNGOSCOPES



ROD LENSES SYSTEM



STAINLESS STEEL BODY



ADAPTERS



RESILIENT SAPPHIRE LENSES



PEEK EYEPIECE

			0°	30°	45°	70°	70°	90°
	Ø	Longitud						
OTOSCOPE	1,9 mm	100 mm	●	●				
	2,7 mm	75 mm	●	●		●		
	4 mm	50 mm	●	●				
RHINOSCOPE	2,7 mm	140 mm	●	●		●		
	2,7 mm	195 mm	●	●		●		
	3 mm	150 mm	●	●	●	●		
	4 mm	185 mm	●	●	●	●		
LARYNGOSCOPE	4 mm	185 mm					●	
	6 mm	158 mm					●	
	10 mm	158 mm					●	●

● Non-Autoclavables ● Autoclavables

- Rod lenses system and stainless steel sealed by laser welding.
- Adapters for connection to type *Storz/Olympus, Wolf* and *ACMI* fiber optic cables.
- Autoclavable endoscopes with lenses made of sapphire of high resistance to scratches.

### OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain  
 Tel: +34 91 804 44 44  
 export@optomic.com  
 www.optomic.com



Specifications are subject to change without notice

