

# DENTAL MICROSCOPE

## OP-Dent 5 Training

SPECIALLY DESIGNED FOR MICROSCOPY TRAINING

- Specific for training and laboratory
- Compact and versatile
- Great mobility in a reduced space



producing quality

The OP-Dent 5 Training has been designed to meet the requirements of microscopy and laboratory training centres. This system allows for the optimization of space, as its arm can be folded to minimum dimensions. Its head comes with integrated LED illumination of high power and long life span.

# DENTAL MICROSCOPE **OP-Dent 5 Training**



Light intensity adjustment



High power LED light



Inclinable binocular (optional)



Instrument tray (Optional)

|                        |  |
|------------------------|--|
| HEAD                   | Stereoscopic head with micrometric focus   |
| BINOCULAR              | 45° inclined (standard), 0° - 240° inclinable and straight (optional)  |
| ILLUMINATION SYSTEM    | LED integrated in head   |
| EYEPIECES              | 10x widefield, 12,5x (optional) and 16x (optional in the inclined binocular)   |
| MAGNIFICATION CHANGER  | 5-steps Galilei  |
| OBJECTIVES             | 250 mm. Optional: 200 mm, 300 mm and 400 mm  |
| HEAD AND ARM           | With brakes and movements using angular contact bearings and high resistance conical cylinder bearings, Soft Move system |
| SUPPORTS               | Table mount  |
| VOLTAGE                | 100/240V 50/60Hz   |
| DIRECTIVE              | Complies with the directive 2017/745 (medical use)   |
| ACCESSORIES (Optional) | Monitor support<br>Instrument Tray<br>Video systems<br>Video-image capture software Digitally Soft                       |



OP-TV8



Compact video microscopy system with built-in CCD high resolution camera.

OP-TV10HD



Compact video microscopy system with built-in CMOS HD high resolution camera.

OP-TV DC



With an exclusive system to adapt it to digital photographic cameras.

## OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain  
Tel: +34 91 804 44 44  
export@optomic.com  
www.optomic.com

