

PORTABLE LED LIGHT SOURCE

OP-LS5

LED TECHNOLOGY MAKES IT EASIER

- Wireless
- Lightweight
- Pure white light



OPTOMIC
producing quality

Powerful LED light source for effective illumination of all kinds of rigid and flexible endoscopes. A convenient alternative for quick and simple inspections.

PORTABLE LED LIGHT SOURCE **OP-LS5**



WOLF TYPE ADAPTER




BATTERY CHARGER AND TRANSFORMER

- Compatible with Storz and Wolf systems.
- LUXEON LED true-white light.
- LED life span: 50.000 hours.
- Ion lithium rechargeable battery.
- Recharge time: 2,5 - 3 hours.

- Autonomy of 3 hours.
- Constant color temperature.
- Weight: 79 gr. (without battery).
- Size: Ø 32 mm x 113 mm.
- Made of aluminium.



RIGID ENDOSCOPE



FLEXIBLE ENDOSCOPE

OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain
Tel: +34 91 804 44 44
export@optomic.com
www.optomic.com



producing quality

Specifications are subject to change without notice

