

# ENT WORKSTATION

## OTOSMART COMPACT

OPTIMIZING SPACE  
IN THE ENT OFFICE

- Compact and modern design
- High capacity shelves and drawers
- Integrated suction and pressure systems



The OTOSMART COMPACT is OPTOMIC's smallest workstation unit. The reduced dimensions makes it the perfect solution for a small practice or as a mobile element between different rooms.



SUCTION AND PRESSURE SYSTEMS



GLOVE BOX SUPPORT AND SPRAY BOTTLE SUPPORT (OPTIONAL)



1L SECRETION CANISTER

## GENERAL CHARACTERISTICS

SUCTION	70L/min - Adjustable and measured by vacuum meter Operating time limit: 60 min. - 1Liter secretion canister
PRESSURE	55L/min - Adjustable power level: 2bars Operating time limit: 30 min.
DESK SURFACE	68 x 50cm - Coloured tempered glass cover included
MAIN DRAWERS	1st and 2nd drawers: 58 x 38 x 6cm (Optional stainless steel tray) 3rd drawer: 58 x 38 x 8cm (Optional stainless steel tray) 4th drawer: 58 x 38 x 13cm
STORAGE SPACE	Bottom storage 39 x 30 x 47 cm with adjustable height shelf
UPPER SHELF	60 x 15 x 1cm - Coloured tempered glass cover included
COLOURS	Blue (other colours available)
MATERIALS	Carbon steel structure / Electrostatic epoxy paint Tempered safety glass
DIMENSIONS WITHOUT ACCESSORIES	121 x 71 x 60cm (87cm Desk height)
WEIGHT WITHOUT ACCESSORIES	133Kg
ELECTRIC OUTPUT	4 Sockets for 230V (120V optional)
MAINS CONNECTION	230V - 50Hz (120V - 60Hz optional)

## OPTIONAL ACCESSORIES

SUPPORTS	1 x Glove box support 1 x Holding up to 3 spray bottles & 4 extra spaces
DRAWER LOCK	The rear key lock blocks all the drawers
CONTROL PEDAL	2 pedals that independently control the activation of suction and pressure

### OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain  
Tel: +34 91 804 44 44 • Fax: +34 91 803 84 14  
export@optomic.com • www.optomic.com

