

ENDOSCOPIC FURNITURE

PROFESSIONAL TROLLEY

YOUR ENDOSCOPIC EQUIPMENT IN THE MINIMUM SPACE

- Adjustable
- Light design
- Easy to handle



Compact and versatile trolley with monitor support, specially designed for your endoscopic equipment using minimum space.
The perfect complement for everyday endoscopic procedures.

ENDOSCOPIC FURNITURE

PROFESSIONAL TROLLEY



Pull-out tray for keyboard and mouse (optional)



Extra monitor support (optional)



Extendable tray with drawer (optional)

- 4 trays; 3 of them are height adjustable in steps of 25 mm.
- Strong shelves: 50 x 60 cm. Maximum load capacity: 10 kg
- Monitor mounting system VESA. 75 & 100.
- Pull-out tray for keyboard and mouse.
- 4 swivel wheels Ø 11,5 cm, two of them with brakes.
- Groove for cables.
- Dimensions: 120 x 71 x 67,5cm cm.
- Weight: 75 kg.
- Colour: light grey.
- Options:
 - Pull-out tray for keyboard and mouse.
 - Monitor articulated arm.
 - Extendable tray with drawer.
 - Holders for 3 rigid endoscopes, videocamera and fiber optic cable.
 - Support for flexible endoscope.

OPTOMIC ESPAÑA S.A.

Madroño 4. E-28770, Colmenar Viejo. Madrid - Spain
Tel: +34 91 804 44 44
export@optomic.com
www.optomic.com

